

Helmer Compartmental Access Refrigerator - iBX080-GX Setup Guide

Location and Placement

NOTICE:


- The Compartmental Access Refrigerator must not be placed in ATEX¹ classified zones as per Directive 99/92/EC ('ATEX 137') and Directive 94/9/EC ('ATEX 95').
 - The Compartmental Access Refrigerator should not be placed in Group 2 medical rooms (ref CEI 64-8 Standard, part 7)
 - The Compartmental Access Refrigerator is classified as IP20 and is not fit for operation outdoors or in environments not protected against atmospheric agents.
- ◆ Has a grounded outlet meeting the electrical requirements stated on the product specification label.
 - ◆ Meets the limits specified for ambient temperature (15°C to 32°C) and relative humidity (80% for temperatures up to 31°C, decreasing linearly to 50% at 40°C).
 - ◆ Is clear of direct sunlight, high temperature sources, heating vents, and air conditioning vents.
 - ◆ Minimum 8" (203 mm) above, and minimum 1" (26 mm).
 - ◆ Install pull handles on cabinet exterior for ease of placement. (optional, if not installed)
 - ◆ Ensure doors are secured and casters are unlocked.
 - ◆ Roll refrigerator into place and lock casters.
 - ◆ Ensure unit is level.
 - ◆ Ensure trays are locked in place within the compartments.
 - ◆ Connect the AC power cord to the receptacle on top of the unit.
 - ◆ Install and secure Access Console (ACX002)
 - ◆ Install and secure Consumable Cart (CCX002) (optional).



(For illustration only; regulatory information and other content shown here may differ from that on the equipment label)

Connect to BloodTrack® Access Console

- 1 Connect the Access Console AC power cord to the blue receptacle on the back of the Access Console.
- 2 Connect the Access Console Ethernet cable to the (LAN2) RJ-45 Ethernet port on the back of the Access Console.
- 3 Route both cables through the holes in the Access Console and the top refrigerator compartment.
- 4 Connect the opposite end of the Access Console power cord to the grey receptacle on the refrigerator.
- 5 Connect the opposite end of the Access Console Ethernet cable to port 6 on the Ethernet switch.



For product registration and complete instructions, refer to current manuals available at

info.helmerinc.com/manuals




Prepare Probe Bottle

- 1 Add glycerin packet to bottle and fill bottle with water.
- 2 Tightly screw cap on probe bottle, place bottle in holder, and insert probe.




Power Unit On

Note: When the refrigerator is first powered on, the Calibration screen will be displayed. The Calibration screen is not displayed on subsequent power-on events.

- 1 Locate the main power switch at the top of the unit and switch it **ON**.
- 3 Switch the i.C³ monitoring/Access Control backup battery switch **ON**. The Start screen is displayed when the i.C³ is powered on.
- 4 If the high temperature alarm sounds, temporarily mute the alarm by touching  Mute button.
- 5 On the Language screen, touch the Language button, then select the preferred language from the drop-down menu.
- 6 Switch Access Console backup battery switch **ON**.

Main Power Switch



Backup Battery Switch

Access Console
Backup Battery Switch



Configure Unit (if necessary)

- ◆ All Helmer refrigerators are preconfigured at Helmer and settings do not need to be changed.
- ◆ Improper adjustment to settings can cause unit failure.
- ◆ If factory default settings need to be adjusted, contact BloodTrackSupport@haemonetics.com or call BloodTrack® Support at 1.877.996.7877.

Tray Operation

- ◆ Trays remain locked at all times unless unlocked by the BloodTrack® system.
- ◆ Trays will be illuminated individually and unlocked when specified via the BloodTrack® kiosk.
- ◆ Under normal operation, trays will lock upon full insertion.
- ◆ A bumper will impede the tray from full extension and removal.
- ◆ Trays have been designed to contain leaks.

