

Model

i.Series® PC2200i Floor Model Platelet Incubator

Specifications

Application, Rating and Electrical Data

Application	Provides the controlled temperature environment required for the storage of platelet products
Agitator Capacity	2 PF96i
Temperature Control Range	20° to 35°C
Electrical Power	115V 50/60Hz 230V 50/60Hz
Fuses	10.0A (230V, quantity 2)
Maximum Current	11.5 FLA / 7.0 FLA
Building Supply Rating	15 amp dedicated ground circuit
Certification / Agency Listing	QPS (Certified to UL and CSA Standards)
Indoor/Outdoor Use	Indoor use only
Application Environment	Non-corrosive, non-flammable, non-explosive, 15°C to 35°C (59°F to 95°F)

Refrigeration System

Refrigeration System	Forced-air circulation
Compressor	1.8 HP (115V) 2.1 HP (230V)
Refrigerant Charge	R134A Non-CFC

Performance Data

Uniformity	+/-1°C
------------	--------

Warranty

Rel.i™ (US/Canada)	5 years compressor, 2 years parts, 1 year labor
Outside US/Canada	Consult your local representative

Controller

Interface	i.Center® LCD viewscreen
Power Switch	On/Off - All Circuit breaker - 230V only
Controller Type	Digital microprocessor
Security	Password protected settings
High / Low Alarms	Fully adjustable
High / Low Alarm Test	Automatic (Peltier)
Door Ajar Alarm	Fully adjustable
Power Failure Alarm	Yes
Condenser Temperature Alarm	Fully adjustable
Event Log	Alarm conditions
AgiTrak System	Yes
Battery Back-up	Yes

Dimensions and Construction

Interior (w x h x d)	36.5 x 35.25 x 18 in 928 x 896 x 458 mm
Exterior (w x h x d)	40.25 x 60 x 28 in 1023 x 1524 x 712 mm
Overall Exterior (w x h x d) (Includes door handle)	40.25 x 60 x 30.25 in 1023 x 1524 x 768 mm
Insulation	1" (26mm)
Exterior / Interior Finish	Bacteria-resistant powder coating
Doors	Dual-pane tempered glass
Key Lock	Yes
Auto-Stop/Start	Yes
Chart Recorder	Yes
Casters	Standard / swivel locking
Net Weight	393 lb (179 kg) 115V 415 lb (189 kg) 230V
Shipping Weight	500 lb (227 kg) 115V 522 lb (237 kg) 230V
Clearance Requirements	Minimum 4" (102mm) across sides and back for ventilation and power cord
Options / Accessories	Remote Alarms, Extended Warranty



