

Storage Trays

Storage Trays from Helmer simplify and organize your storage needs. Trays are available for all Helmer refrigerators and freezers including Undercounter, Upright, and Pass-Thru models.

Storage tray packages come with dividers and rails so that each storage area can be customized. The dividers are removable and adjustable to accommodate changing needs.

Shown with RT1000-P in roll out basket

Packages

RT1000-P Storage Tray Half-Size
For use with all models
Includes 1 rail and 5 dividers
10" x 15 1/4" x 2 5/16"
(254mm x 388mm x 59mm)

RT0511-P Storage Tray Full-size
For use with U/C and 11 cf models
Includes 2 rails and 6 dividers
15 1/4" x 16" x 2 5/16"
(388mm x 407mm x 59mm)

RT2000-P Storage Tray Long Rectangle
For use with 25 and 56 cf models
Includes 5 dividers
5 1/2" x 26" x 2 5/16"
(140mm x 661mm x 59mm)

RT2045-P Storage Tray Full-size
For use with 20 and 45 cf models
Includes 3 rails and 12 dividers
20" x 20" x 2 5/16"
(508mm x 508mm x 59mm)

RT2556-P Storage Tray Full-Size
For use with 25 and 56 cf models
Includes 3 rails and 12 dividers
20" x 26" x 2 5/16"
(508mm x 661mm x 59mm)


Individual trays also sold separately



RT-1000-P




RT-0511-P



RT-2000-P



RT-2045-P



RT-2556-P

