

Quality in motion.

## TechSpin™ from Helmer

### Cost Effective, Flat Separation



The EBA 270 benchtop centrifuge offers a 90° swing-out in a small footprint for flat separation and optimal results.

- Accepts common blood and urine tubes up to 15ml
- 6 easily removable / interchangeable tube carriers
- Programmable in RPM and minutes
- Smooth plastic housing, easy to clean
- Stainless steel, seamless chamber



### Package Specifications

Max RPM / RCF	4,000 / 2,254
Max Package Volume	90ml
Spin Capacity	6 13 x 75/100 6 16 x 75/100
Electrical	115V / 50-60Hz
Dimensions	13" x 9.5" x 15.5"
Weight	29 lb / 13 kg

### TechSpin™ Package

Part No. HET2300-01

Includes:

- EBA 270 centrifuge
- 2315 6-place swing-out rotor
- 2331 insert, 15 ml (6 each)
- 2333 insert, 10 ml (6 each)



Helmer® and TechSpin™ are trademarks of Helmer, Inc. in the USA  
© 2011 Helmer, Inc. 380239-1/A 4/11



14395 Bergen Boulevard • Noblesville, IN 46060  
Toll Free (US & Can): 800.743.5637  
Ph: 317.773.9073 • Fax: 317.773.9082  
sales@helmerinc.com • www.helmerinc.com