

Quality in motion.

Microcentrifugation from Helmer

Options for Managing Your Microcentrifuge Needs

Choose from our 3 benchtop models that range in capacity from 12 to 48 microtubes. These quality centrifuges offer many features for the critical needs of your lab.

- Powerful performance
- Bio-containment options
- Multiple rotor and tube options
- Reliable, quiet operation



Mikro 120

- Compact footprint
- 24 place rotor
- Quick and easy to use
- Impulse key for short spins



Mikro 200 / 200R

- Available in non-refrigerated and refrigerated versions
- Multiple microliter rotors available:
 - 24 or 30 place rotors
 - PCR strip rotor
 - Cryo tube rotor
- Bio-containment or phenol-resistant lids
- Fast Cool function (200R)



Mikro 220 / 220R

- Available in non-refrigerated and refrigerated versions
- Multiple microliter rotors available:
 - 12, 24, or 30 place rotors
 - High capacity 48 place rotor
- Mikro 220R offers extended range of applications with additional rotors for blood tubes, conical tubes, and more

Microcentrifuges and Rotors

Mikro 120				
Part #	# Places	RPM	RCF	Lid
1242	24-place	14,000	18,407	metal
Mikro 200/200R				
Part #	# Places	RPM	RCF	Lid
2424-B	24-place	15,000	21,382	plastic optional: autoclavable bio-containment or autoclavable bio-containment phenol resistant lids
2427	30-place, 2 row	15,000	20,627	plastic optional: autoclavable bio-containment or autoclavable bio-containment phenol resistant lids
2430-A	30-place	15,000	20,376	plastic optional: autoclavable bio-containment or autoclavable bio-containment phenol resistant lids
2426-A	20-place, cryo	15,000	21,130	plastic optional: autoclavable bio-containment or autoclavable bio-containment phenol resistant lids
2418	4-place, PCR strip	15,000	14,338	optional: plastic lid
Mikro 220/220R				
Part #	# Places	RPM	RCF	Lid
2218	12-place	18,000	25,718	plastic
1195-A	24-place	18,000	31,514	metal bio-containment lid
1189-A	30-place	14,000	21,255	metal bio-containment lid
1158-L	48-place, 2 row	14,000	21,255 outer 18,845 inner	metal bio-containment lid
2219	20-place, cryo	14,000	18,407	plastic optional: autoclavable bio-containment or autoclavable bio-containment phenol resistant lids
1154-L	24-place, swing-out	13,000	18,516	bio-containment lid
1161	6-place, drum rotor	13,000	14,171	metal lid
1160	6-place, PCR strip	14,000	18,845	optional: plastic lid
1023	24-place, hematocrit	15,000	21,382	evaluation disk lid
1015*	12-place, fixed angle for blood tubes	6,000	4,146	none
1016*	6-place, fixed angle for 15ml / 50ml conical tubes	6,000	4,025	none
2226*	12-place, swing out for blood tubes	5,000	2,963	none
1020*	8-place, swing out for blood tubes	5,000	2,879	none

* Available on 220R only

