































































	Tubes 10/12 x 75mm	Conicals 15ml 50ml		Tubes and Bottles	Microtubes 1.5- 2.0ml	Shallow	Deep Well	Blood Collection Tubes	Hematocrit Capillaries	Blood Bags	Special Tubes	Cyto Applications
Small Centrifuges												
EBA 20	8 each 	4 each 						8 each 				
EBA 20S	8 each 	4 each 						8 each 				
EBA 21	12 each 	6 each 	3 each 	6 x 50ml 	30 each 			12 each 	24 caps 			
Benchtop Centrifuges												
Rotofix 32A	68 each 	32 each 	8 each 	4 x 100ml 	36 each 			28 each 			4 each 	6 x 8ml 
Universal 320 / 320R	68 each 	32 each 	8 each 	4 x 100ml 	48 each 	10 each 	2 each 	32 each 	24 caps 		4 each 	6 x 8ml 
Rotina 38 / 38R	102 each 	28 each 	10 each 	4 x 250ml 	160 each 	10 each 	2 each 	48 each 			4 each 	8 x 8ml 
Rotina 420 / 420R	136 each 	52 each 	20 each 	4 x 600ml 	160 each 	16 each 	4 each 	84 each 			4 each 	8 x 8ml 
Rotanta 460 / 460R 460RS	160 each 	56 each 	24 each 	4 x 750ml 	320 each 	24 each 	4 each 	96 each 		4 each* 	6 each 	12 x 8ml 

*Rotanta 460RS accommodates 4 blood bags (450ml)