

# NEW! Undercounter Freezers

Offer Secure Cold Storage Solutions



## i.Series® and Horizon Series™

- Laboratory and Plasma models available
- Programmable defrost times
- Undercounter, free-standing and stackable
- Exclusive Sure-Seal door design
  - Self-closing
  - Magnetic closure eliminates broken door latches
  - Field reversible
  - Key lock
- Bacteria-resistant powder coating
- Customize with: swivel casters, roll-out baskets, and stainless steel interior/exterior

**Now Accepting Orders!**

## Specifications

Dimensions (wxhxd)	24" x 33.5" x 26.75"
Interior Volume	5 cu ft
Operating Temperature Range	-10°C to -30°C



[www.helmerinc.com](http://www.helmerinc.com)

800-743-5637