

Storage Trays

Molded trays with flexible dividers

Storage Trays simplify and organize storage of medication for standardization, security, and control. Trays are available for all Helmer refrigerators and freezers and can be used in stainless steel drawers, roll-out baskets, and shelves.


Storage tray packages come with dividers and rails so that each storage area can be customized for maximum efficiency. The dividers are removable and adjustable to accommodate changing needs.




Packages

Standard Tray & Divider Set-up

#400891-1 U/C & 11 cf
Includes 2 rails and 6 dividers
15 1/4" x 16" x 2 5/16"
(388mm x 407mm x 59mm)



#400892-1 20 & 45 cf
Includes 3 rails and 12 dividers
20" x 20" x 2 5/16" (508mm x 508mm x 59mm)



#400893-1 25 & 56 cf
Includes 3 rails and 12 dividers
20" x 26" x 2 5/16" (508mm x 661mm x 59mm)


Small Tray & Divider Set-up

#400894-1 All models
Includes 1 rail and 5 dividers
10" x 15 1/4" x 2 5/16"
(254mm x 388mm x 59mm)



Long Tray & Divider Set-up

#400895-1 25 & 56 cf
Includes 5 dividers
5 1/2" x 26" x 2 5/16"
(140mm x 661mm x 59mm)



Trays and dividers also sold separately

www.helmerinc.com | 800.743.5637

380223-1/D 11/15 © 2015 Helmer, Inc. All rights reserved.