



CENTRIFUGES

EBA 280

SEROLOGY / STAT / PEDIATRIC

Separating Ourselves 90° From the Competition

High-performance STAT centrifuge with 90° horizontal rotor for flat separation in all standard blood, urine and pediatric tubes. Achieve consistent platelet counts in minutes with the EBA 280 S model.

MAX. CAPACITY	6 x 50 mL
MAX. RPM/RCF	6,000/5,071

Distributed by:



QUIET OPERATION

Quiet as a normal conversation. Will not increase noise levels or contribute to distractions.



FIVE YEAR WARRANTY

All Hettich centrifuges are protected by our 5-year warranty and full technical support.



FUNCTIONAL DESIGN

Easy to lift lid with front-facing control panel for effortless interaction and programming.



FULLY CERTIFIED

Engineered and manufactured in accordance with all relevant safety and production standards.



ENHANCED SAFETY

Self-closing, auto-locking lid with steel latches. All-metal chamber provides reinforced security.



MADE IN GERMANY

German engineered and manufactured to be quiet, reliable and safe since 1904.

DIMS. (HxWxD): 9.5" x 12.8" x 15.75"

→ hettweb.com/eba280

HETTICH EBA 280 APPLICATION PACKAGES



EBA 280 BLOOD TUBE PACKAGE - 6

	Tube Qty.	Volume (mL)	Size (mm)	280	280 S	280	280 S	PKG. NO. (280)	PKG. NO. (280 S)
				RPM		RCF			
	6	1.6 - 5	13 x 75	≤ 4,700	≤ 6,000	≤ 2,618	≤ 4,266	280BLOOD6	280SBLOOD6
	6	4 - 7	13 x 100	≤ 4,700	≤ 6,000	≤ 3,112	≤ 5,071		
	6	4 - 7	16 x 75	≤ 4,700	≤ 6,000	≤ 2,865	≤ 4,669		
	6	8.5 - 10	16 x 100	≤ 4,700	≤ 6,000	≤ 3,112	≤ 5,071		
Includes: (1) Centrifuge, (1) Horizontal Rotor, (6) Buckets, (6) Adapters 13 x 75 mm, (6) Adapters 16 x 75 mm									

EBA 280 BLOOD TUBE PACKAGE - 8

	Tube Qty.	Volume (mL)	Size (mm)	RPM		RCF		PKG. NO.
				280	280 S	280	280 S	
	8	1.6 - 5	13 x 75	≤ 5,000		≤ 2,991		280BLOOD8
	8	4 - 7	16 x 75	≤ 5,000		≤ 2,991		
Includes: (1) Centrifuge, (1) Horizontal Rotor, (8) Buckets 13 x 75 mm, (8) Buckets 16 x 75 mm								

EBA 280 BLOOD TUBE PACKAGE - 12

	Tube Qty.	Volume (mL)	Size (mm)	RPM		RCF		PKG. NO.
				280	280 S	280	280 S	
	12	1.6 - 5	13 x 75	≤ 6,000		≤ 3,300		280BLOOD12
	12	4 - 7	13 x 100	≤ 6,000		≤ 4,146		
	12	4 - 10	16 x 75/100	≤ 6,000		≤ 4,146		
Includes: (1) Centrifuge, (1) Angle Rotor, (12) Adapters 13 x 75/100 mm								

EBA 280 PEDIATRIC TUBE PACKAGE - 6

	Tube Qty.	Volume (mL)	Size (mm)	280	280 S	280	280 S	PKG. NO. 280	PKG. NO. 280 S
				RPM		RCF			
	6	0.5	10.7 x 36	≤ 4,700	≤ 6,000	≤ 1,877	≤ 3,059	280BLOOD6-PED	280SBLOOD6-PED
	Includes: (1) Centrifuge, (1) Horizontal Rotor, (6) Buckets, (6) Adapters								

EBA 280 SEROLOGY PACKAGE - 12

	Tube Qty.	Volume (mL)	Size (mm)	RPM		RCF		PKG. NO.
				280	280 S	280	280 S	
	12	1.6 - 5	10/12/13 x 75	≤ 5,000		≤ 2,963		280SEROLOGY
	Includes: (1) Centrifuge, (1) Angle Rotor, (12) Buckets							

EBA 280 SEROLOGY DECANTING PACKAGE - 12

	Tube Qty.	Volume (mL)	Size (mm)	RPM		RCF		PKG. NO.
				280	280 S	280	280 S	
	12	5 - 7	12 x 75/100	≤ 5,000		≤ 2,879		280SEROLOGY-D
	12	5 - 7	10 x 75	≤ 5,000		≤ 2,879		
Includes: (1) Centrifuge, (1) Angle Rotor								

EBA 280 URINALYSIS PACKAGE - 6

	Tube Qty.	Volume (mL)	Size (mm)	RPM		RCF		PKG. NO.
				280	280 S	280	280 S	
	6	12 KOVA	17 x 102	≤ 4,700		≤ 3,112		280URINALYSIS6
	Includes: (1) Centrifuge, (1) Horizontal Rotor, (6) Buckets							

EBA 280 CUSTOM PACKAGE

	Supports over 64 application configurations, please inquire for more information								
Tube types include: Blood, Urine, Pediatric, Standard and Conical up to 50 mL									



hettweb.com/eba280

DISTRIBUTED BY:	
	Helmer Scientific 14400 Bergen Boulevard, Noblesville, IN 46060 USA Toll Free: (800) 743-5637 Email: sales@helmerinc.com www.helmerinc.com



Discover more solutions online at hettweb.com **Toll Free:** 1(866) 370-4388 **Email:** info@hettweb.com
 © 2016 Hettich Instruments, LP. All rights reserved. Lit. No. PS-EBA-280-E Rev A 01-2016.

