


Drawer Dividers



Drawer Dividers from Helmer simplify and organize your storage needs. They are excellent for storing blood bags and for use with automated medication dispensing systems.

Drawer Dividers are designed exclusively for Helmer and fit conveniently in our stainless steel drawers. They are available in angled or upright configurations. The heavy-duty plastic divider parts are removable to create custom storage spaces.



Upright Drawer Dividers

	Compartments	Compartment Dimensions (w x d) (in/mm)
5/11 cf (142/326 liter)	28	4.1 x 2.3 (105 x 59)
20/45 cf (572/1271 liter)	40	4.2 x 2.4 (107 x 61)
25/26.5/56/58 cf (714/750/1586/1642 liter)	55	4.2 x 2.3 (107 x 59)



Angled Drawer Dividers

	Compartments
5/11 cf (142/326 liter)	24
20/45 cf (572/1271 liter)	30
25/26.5/56/58 cf (714/750/1586/1642 liter)	40