

Locking Drawer

For Refrigerators

Locking drawers offer isolated storage for your medications, vaccines, blood products, and more.¹ The ventilated cover is secured by a key lock adding an extra layer of protection while maintaining excellent temperature uniformity.

Locking drawers can be used in undercounter and upright refrigerator models in both i.Series® and Horizon Series™ lines.



Features

- Available for undercounter and upright models (4-56 cu ft)
- Key lock
- Solid stainless steel pull out drawer with ventilated cover
- Scratch-Guard edge protectors keep glass doors from being scratched by drawers



Specifications

	Depth	Width	Height	Drawers*
4 and 5 cf	16.75" / 426mm	17" / 432mm	4" / 102mm	2
11 cf	16.75" / 426mm	17" / 432mm	4" / 102mm	5
20 cf	22" / 559mm	21" / 534mm	4" / 102mm	1*
45 cf	22" / 559mm	21" / 534mm	4" / 102mm	2*
25 cf	22" / 559mm	27" / 636mm	4" / 102mm	1*
56 cf	22" / 559mm	27" / 686mm	4" / 102mm	2*

* Must be placed in bottom position(s)