

Pass-Thru Refrigerator Accessories

Finish Trim Kit


The Finish Trim Kit provides a perimeter frame around the Pass-Thru Refrigerator. It adheres to the walls surrounding the refrigerator and provides a superior, finished look that is easy to clean.

Finish Trim Kits are bacteria-resistant powder coated to match the Pass-Thru Refrigerator. No hardware is required for installation. Kits are available for single and double door models.

Skirt Kit

The Skirt mounts to the bottom of the refrigerator on the clean room side sealing the space between the refrigerator and floor.

The Skirt is bacteria-resistant powder coated to match the Pass-Thru Refrigerator. Hardware is provided for installation. Skirts are available for single and double door models.



Finish Trim Kit

Includes 2 side rails, 1 top rail, and VHB mounting tape. To enclose both sides of the refrigerator, 2 Finish Trim Kits are needed.

Single Door

Side Rails 3" x 81.75" / Top Rail 3" x 35"

Double Door

Side Rails 3" x 81.75" / Top Rail 3" x 65"

Multiple Unit Kit - 2 Single Doors

Side Rails 3" x 81.75" / Top Rail 3" x Adjusts to fit two single door refrigerators placed side by side (total width 59")

Multiple Unit Kit - 1 Single and 1 Double Door

Side Rails 3" x 81.75" / Top Rail 3" x Adjusts to fit one single door refrigerator and one double door refrigerator placed side by side (total width 88.75")

Skirt Kit

Includes skirt and installation hardware. Available for clean room side only.

Finish Trim Kit and Skirt Kit sold separately