


# Temperature Chart Recorders

## TCR1 / TCR2 / TCR3

The three TCR models provide accurate and reliable temperature recording. These durable standalone units are versatile and easy to use. Each TCR model continuously records the temperature within a specified operating range. In addition, they offer complete flexibility as each unit may be programmed to run at any of the three available ranges.




### Features

- 4" (102mm), 7-day circular standalone chart recorder with inkless, pressure sensitive chart paper
- Countertop, wall, or under cabinet placement
  - Stand with tilt adjustment knob included for countertop or under cabinet use
  - Mounting holes for wall placement
- Battery backup with visual low battery indicator
- 15 ft (5m) cable with RTD probe
- Durable, galvanized steel housing with bacteria-resistant powder coating
- User programmable to operate at any of the three operating ranges
- Included with the TCR: stylus, battery, 15 ft (5m) probe, 1 package of 52 charts, 7 ft (2.5m) power cord, and Operation Manual
- Certified to applicable UL and CSA standards by a NRTL


### Mounting Options



COUNTERTOP



WALL



UNDER CABINET

### Temperature Ranges

Model	Application	Operating Range
TCR1	Incubator/Ambient	0 to 35° C
TCR2	Refrigerator	-5 to 20° C
TCR3	Freezer	-50 to 0° C

### Chart Recorder Paper

4" (102mm) diameter, package of 52

TCR1	Part No. 220273	0 to +35° C
TCR2	Part No. 220366	-5 to +20° C
TCR3	Part No. 220419	-50 to 0° C

### Specifications

Dimensions w x h x d	Weight	Electrical
8" x 7" x 4.5"	5.1 lbs / 2.3 kg	115V 50/60Hz
204mm x 178mm x 115mm		230V 50/60 Hz

Dimensions with stand: 9.5" w x 8.25" h x 6.75" d (242mm x 210mm x 172mm). Overall exterior dimensions for the chart recorder with the stand, adjusted so chart recorder is tilted about 30° and including items that protrude from the main unit.