

HETTICH CLINICAL TABLETOP CENTRIFUGES

EBA SERIES REBATE PROMOTION



SAVE

FEATURES:

Purchase a Hettich **EBA 20, EBA 20S, EBA 270 or EBA 21** Centrifuge Package for Blood or Urinalysis and save. Offer valid for a limited time. See your local rep for details.

- \$100 Customer Rebate
- Free Shipping and Handling
- FDA registered manufacturer
- Run-Silent™ Operation (as quiet as normal conversation)
- Safety-Lock™ Lid
- Five (5) Year Warranty

Hettich
CENTRIFUGES



EBA 20

BLOOD/ URINALYSIS

- Built-in, multi-functional, 8-place angled rotor
- Max. Capacity 8 tubes
- Max. RPM/ RCF: 6,000/ 3,461

EBA 20S

STAT/ COAG./ PPP

- 3-minute spin for platelet counts less than 10,000/ μ L
- Built-in, multi-functional, 8-place angled rotor
- Max. Capacity 8 tubes
- Max. RPM/ RCF: 8,000/ 6,153
- Maintenance free/ brushless motor

EBA 21

BLOOD/ SEROLOGY/ URINALYSIS

- 90° Swing-out & Fixed angle rotors available
- Max. Capacity: 12 blood tubes
- Max. RPM/RCF: 18,000/ 25,718
- Maintenance free/ brushless motor

EBA 270

GEL BLOOD TUBES

- Built-in 90° swing-out rotor
- Max. Capacity (6) 16x100mm
- Max. RPM.RCF: 4,000/2,254
- Ideal for the physicians office
- Maintenance free/ brushless motor




Download the rebate form by visiting:
www.hettweb.com/eba-rebate-form.pdf


**HETTICH PARTNER
PROMOTION**

PROMOTION VALID THROUGH PARTICIPATING DISTRIBUTION PARTNER ONLY.


EBA 20 BLOOD TUBE PACKAGE

	Tube Qty	Volume (mL)	Size (mm)	RPM	G-FORCE	PACKAGE CONTENTS
8	1.6 - 7	13 x 75/ 100	6,000	3,461	2,697	(1) Centrifuge (1) Fixed Angle Built-in Rotor (8) Adapters for 13 x 75/ 100 mm (8) Adapters for 16 x 100 mm
	8.5 - 10	16 x 100				
Pkg No.:				20BLOOD		


EBA 20 URINALYSIS PACKAGE

	Tube Qty	Volume (mL)	Size (mm)	RPM	G-FORCE	PACKAGE CONTENTS
4	KOVA (12)	17 x 102	6,000	3,461	(1) Centrifuge (1) Fixed Angle Built-in Rotor	
	10 Conical	--				
Pkg No.:				20URINALYSIS		


EBA 20S BLOOD TUBE PACKAGE

	Tube Qty	Volume (mL)	Size (mm)	RPM	G-FORCE	PACKAGE CONTENTS
8	1.6 - 7	13 x 75/ 100	8,000	6,153	4,794	(1) Centrifuge (1) Fixed Angle Built-in Rotor (8) Adapters for 13 x 75 mm (8) Adapters for 13 x 100 mm
	8.5 - 10	16 x 100				
Pkg No.:				20SBLOOD		


EBA 21 BLOOD TUBE PACKAGE - 8

	Tube Qty	Volume (mL)	Size (mm)	RPM	G-FORCE	PACKAGE CONTENTS
8	1.6 - 7	13/ 16 x 75	5,000	2,879	(1) Centrifuge (1) 90° Swing-out Rotor (8) Buckets for 13 x 75 mm (8) Buckets for 16 x 75 mm	
	Pkg No.:					21BLOOD8


EBA 21 BLOOD TUBE PACKAGE - 12

	Tube Qty	Volume (mL)	Size (mm)	RPM	G-FORCE	PACKAGE CONTENTS
12	1.6 - 7	13 x 75/ 100	6,000	4,146	3,300/ 4,146	(1) Centrifuge (1) Fixed Angle Rotor (12) Adapters for 13 x 75 mm (12) Adapters for 13 x 100 mm
	8.5 - 10	16 x 100				
Pkg No.:				21BLOOD12		

EBA 21 URINALYSIS PACKAGE

	Tube Qty	Volume (mL)	Size (mm)	RPM	G-FORCE	PACKAGE CONTENTS
6	KOVA (12)	17 x 102	6,000	4,416	(1) Centrifuge (1) Fixed Angle Rotor	
	Pkg No.:					21URINALYSIS

EBA 270 BLOOD TUBE PACKAGE

	Tube Qty	Volume (mL)	Size (mm)	RPM	G-FORCE	PACKAGE CONTENTS
6	1.6 - 7	13 x 75/ 100	4,000	1,807/ 2,254	(1) Centrifuge (1) Swing-out Rotor (6) Adapters for 13 x 75 mm (6) Adapters for 13 x 100 mm	
	8.5 - 10	16 x 75/ 100				
Pkg No.:				270BLOOD		

→ Download the rebate form by visiting: www.hettweb.com/eba-rebate-form.pdf

DISTRIBUTED BY:



Helmer Scientific

14400 Bergen Boulevard, Noblesville, IN 46060

Toll Free: 800-743-5637 or +1-317-773-9073

Email: service@helmerinc.com



Hettich Promotion is valid only through a participating Hettich Distribution Partner. Standard ground shipping is free of charge (excluding holidays and weekends). The Hettich Five (5) Year Warranty adds an additional 3-year powertrain warranty to the Hettich standard Two (2) Year Warranty; restrictions apply, inquire for more details. Offer is valid in the U.S. only, January 01, 2014 to December 31, 2014. No substitutions.

HETTICH NORTH AMERICA

Toll Free: 1-866-370-4388

Email: info@hettweb.com

→ Discover more solutions online at: www.hettweb.com/centrifuges

Copyright © 2012 Hettich Instruments, LP. All rights reserved. For the most up to date information on Hettich please visit hettweb.com.
Printed in the USA using soy based inks on FSC certified recycled paper (25% pre-consumer waste & 25% post consumer waste).

PS-EBA-REBATE-E RevA 01-2014