

Door Handles and Adapters for Pyxis® Lock

Special accessories are available for facilities that install the Pyxis Remote Manager and Smart Remote Manager on Helmer upright refrigerators and undercounter refrigerators.

Pyxis Compatible Door Handles (Upright Models)

The door handle features a cutout and spacer which enable the Pyxis lock to wrap around the front and side of the door.


- Available for single door upright refrigerators 11, 20, 25, and 26.5 cu ft (326, 572, 714 and 750 L)
- Compatible with glass and solid doors
- Factory or field installation

Pyxis Lock Adapter Kit (Undercounter Models)

The adapter kit includes an extension bar and alternate handle that mount to the door of the unit, allowing the Pyxis lock to be properly installed.

- Adapter kits available for undercounter refrigerators 4 and 5 cu ft (113 and 142 L)
- Compatible with glass and solid doors; powder coated and stainless steel exteriors
- Factory installed or field installation kits available

Cut-out of Door Handle for Upright Refrigerators



Upright Pyxis®
Compatible
Door Handle



Undercounter Adapter
(Shown with Pyxis® Lock)



Description

Pyxis Lock Adapter Kit, Factory installed, 4 and 5 cf (113 and 142 L) (i.Series® & Horizon Series™)

Pyxis Lock Adapter Kit, Factory installed, 4 and 5 cf (113 and 142 L) (Scientific Series™)

Pyxis Lock Compatible Door Handle, Factory installed, 11 cf (326 L)

Pyxis Lock Compatible Door Handle, Factory installed, 20, 25 and 26.5 cf (572, 714 and 750 L)

Pyxis Lock Adapter Kit, Field installation, 4 and 5 cf (113 and 142 L) (i.Series® & Horizon Series™)

Pyxis Lock Adapter Kit, Field installation, 4 and 5 cf (113 and 142 L) (Scientific Series™)

Pyxis Lock Compatible Door Handle, Field installation, 11 cf (326 L)

Pyxis Lock Compatible Door Handle, Field installation, 20, 25 and 26.5 cf (572, 714 and 750 L)

Note: Adapter kits and compatible door handles do not include the medication dispensing lock.

