



CENTRIFUGES

EBA 200

BLOOD / URINALYSIS / PEDIATRIC

Press 'Start' in Case of Emergency

High-performance with built-in 8-place rotor, compact frame and all-metal containment for added safety. Ideal for blood, urine and pediatric tubes. STAT, PPP and coag. in minutes with the EBA 200 S model.

MAX. CAPACITY	8 x 15 mL
MAX. RPM/RCF	8,000/6,153

Distributed by:



QUIET OPERATION

Quiet as a normal conversation. Will not increase noise levels or contribute to distractions.



FIVE YEAR WARRANTY

All Hettich centrifuges are protected by our 5-year warranty and full technical support.



FUNCTIONAL DESIGN

Easy to lift lid with front-facing control panel for effortless interaction and programming.



FULLY CERTIFIED

Engineered and manufactured in accordance with all relevant safety and production standards.



ENHANCED SAFETY

Self-closing, auto-locking lid with steel latches. All-metal chamber provides reinforced security.



MADE IN GERMANY

German engineered and manufactured to be quiet, reliable and safe since 1904.

DIMS. (HxWxD): 9.0" x 10.3" x 13.9"

→ hettweb.com/eba200

HETTICH EBA 200 APPLICATION PACKAGES



EBA 200 BLOOD TUBE PACKAGE - 8

Tube Qty.	Volume (mL)	Size (mm)	200	200 S	200	200 S
			RPM		RCF	
8	1.6 - 5	13 x 75	≤ 6,000	≤ 8,000	≤ 2,697	≤ 4,794
8	4 - 7	13 x 100	≤ 6,000	≤ 8,000	≤ 3,461	≤ 6,153
8	8.5 - 10	16 x 100	≤ 6,000	≤ 8,000	≤ 3,461	≤ 6,153
4	8	16 x 125	≤ 6,000	≤ 8,000	≤ 3,461	≤ 6,153

Includes: (1) Centrifuge, (1) Built-in Angle Rotor, (8) Adapters 13 x 75 mm, (8) Adapters 13 x 100 mm

PKG. NO. (EBA 200)	PKG. NO. (EBA 200 S)
200BLOOD8	200SBLOOD8

EBA 200 PEDIATRIC TUBE PACKAGE - 8

Tube Qty.	Volume (mL)	Size (mm)	200	200S	200	200 S
			RPM		RCF	
8	0.5	10.7 x 36	≤ 6,000	≤ 8,000	≤ 2,214	≤ 3,935

Includes: (1) Centrifuge, (1) Built-in Angle Rotor, (8) Adapters 0.5 mL

PKG. NO. (EBA 200)	PKG. NO. (EBA 200 S)
200BLOOD8-PED	200SBLOOD8-PED

EBA 200 URINALYSIS PACKAGE - 4

Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF
4	12 KOVA	17 x 102	≤ 6,000	≤ 3,461
4	15 Conical	17 x 120	≤ 6,000	≤ 3,461

Includes: (1) Centrifuge, (1) Built-in, Angle Rotor

PKG. NO.
200URINALYSIS4

EBA 200 CUSTOM PACKAGE

	Supports over 19 application configurations, please inquire for more information
	Tube types include: Standard, Blood, Pediatric, Urine, and Conical up to 15 mL



hettweb.com/eba200

DISTRIBUTED BY:



Helmer Scientific
 14400 Bergen Boulevard, Noblesville, IN 46060 USA
Toll Free: (800) 743-5637 | **Email:** sales@helmerinc.com
www.helmerinc.com



Discover more solutions on-line at hettweb.com **Toll Free:** 1(866) 370-4388 **Email:** info@hettweb.com
 © 2016 Hettich Instruments, LP. All rights reserved. Lit. No. PS-EBA-200-E Rev A 01-2016.



LAB TECHNOLOGY