



CENTRIFUGES

HAEMATOKRIT 200

HEMATOCRIT CAPILLARIES

Seven Minutes of Separation

Spin up to 24 capillaries per run and achieve best results within seven minutes. Rotor lid doubles as an evaluation disc allowing you to determine hematocrit values without having to remove the capillaries and disturb the sedimentation. Also supports a 20-place rotor for buffy coat analysis.

MAX. CAPACITY	24 Capillaries
MAX. RPM/RCF	13,000/16,060

Distributed by:



QUIET OPERATION

Quiet as a normal conversation. Will not increase noise levels or contribute to distractions.



FIVE YEAR WARRANTY

All Hettich centrifuges are protected by our 5-year warranty and full technical support.



FUNCTIONAL DESIGN

Easy to lift lid with front-facing control panel for effortless interaction and programming.



FULLY CERTIFIED

Engineered and manufactured in accordance with all relevant safety and production standards.



ENHANCED SAFETY

Self-closing, auto-locking lid with steel latches. All-metal chamber provides reinforced security.



MADE IN GERMANY

German engineered and manufactured to be quiet, reliable and safe since 1904.

DIMS. (HxWxD): 9.0" x 10.3" x 13.9"

hettweb.com/haematokrit200


HETTICH HAEMATOKRIT 200 APPLICATION PACKAGES



HAEMATOKRIT 200 HEMATOLOGY PACKAGE - 24

	Capillary Qty.	Volume (mL)	Size (mm)	RPM	RCF	PKG. NO.
	24	Capillaries	--	≤ 13,000	≤ 16,060	200HEMA24
Includes: (1) Centrifuge, (1) Hematocrit Rotor, (1) Lid, Evaluation disc, (1) Sealing Putty						

HAEMATOKRIT 200 HEMATOLOGY BUFFY COAT PACKAGE - 20

	Capillary Qty.	Volume (mL)	Size (mm)	RPM	RCF	PKG. NO.
	20	Capillaries	--	≤ 13,000	≤ 16,060	200HEMA20
Includes: (1) Centrifuge, (1) Hematocrit Disc Rotor, (1) Lid						



hettweb.com/haematokrit200

DISTRIBUTED BY:



Helmer Scientific
 14400 Bergen Boulevard, Noblesville, IN 46060 USA
Toll Free: (800) 743-5637 | **Email:** sales@helmerinc.com
www.helmerinc.com



Discover more solutions online at hettweb.com **Toll Free:** 1(866) 370-4388 **Email:** info@hettweb.com
 © 2016 Hettich Instruments, LP. All rights reserved. Lit. No. PS-HAEMATOKRIT-200-E Rev A 01-2016.



LAB TECHNOLOGY