



CENTRIFUGES

# MIKRO 185

MICROLITER / SPIN COLUMNS / PEDIATRIC

## Another Check in the Spin Column

Offers application versatility by enabling users to run up to 24 microliter tubes, 12 pediatric tubes, or 18 spin column kits in a single, ultra-compact frame. High-performance with all-metal containment for added safety.

<b>MAX. CAPACITY</b>	24 x 1.5/2.0 mL
<b>MAX. RPM/RCF</b>	13,300/17,008

Distributed by:



### QUIET OPERATION

Quiet as a normal conversation. Will not increase noise levels or contribute to distractions.



### FIVE YEAR WARRANTY

All Hettich centrifuges are protected by our 5-year warranty and full technical support.



### FUNCTIONAL DESIGN

Easy to lift lid with front-facing control panel for effortless interaction and programming.



### FULLY CERTIFIED

Engineered and manufactured in accordance with all relevant safety and production standards.



### ENHANCED SAFETY

Self-closing, auto-locking lid with steel latches. All-metal chamber provides reinforced security.



### MADE IN GERMANY

German engineered and manufactured to be quiet, reliable and safe since 1904.


**DIMS. (HxWxD):** 9.0" x 10.3" x 13.9"

→ [hettweb.com/mikro185](http://hettweb.com/mikro185)

# HETTICH MIKRO 185 APPLICATION PACKAGES




## MIKRO 185 MICROLITER TUBE PACKAGE - 12

	Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF	PKG. NO.
	12	1.5 - 2.0	11 x 38	≤ 13,300	≤ 14,041	185MICRO12
	<b>Includes:</b> (1) Centrifuge, (1) Fixed-angle Rotor, (12) Adapters 1.5 mL					

## MIKRO 185 MICROLITER TUBE PACKAGE - 24 - BIO-CONTAINMENT

	Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF	PKG. NO.
	24	1.5 - 2.0	11 x 38	≤ 13,300	≤ 17,008	185MICRO24
	<b>Includes:</b> (1) Centrifuge, (1) Fixed-angle Rotor, (24) Adapters 1.5 mL, (1) Bio-containment Lid					

## MIKRO 185 PEDIATRIC TUBE PACKAGE - 12


	Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF	PKG. NO.
	12	0.5	10.7 x 36	≤ 13,300	≤ 16,019	185MICRO12-PED
	12	0.8	8 x 45	≤ 13,300	≤ 17,008	
	<b>Includes:</b> (1) Centrifuge, (1) Fixed-angle Rotor, (12) Adapters 0.5 mL, (12) Adapters 0.8 mL					

## MIKRO 185 SPIN COLUMN PACKAGE - 18

	Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF	PKG. NO.
	18	1.5 - 2.0	11 x 38	≤ 13,300	≤ 15,030	185MICROSC
	<b>Includes:</b> (1) Centrifuge, (1) Fixed-angle Rotor, (18) Adapters 1.5 mL, (1) Bio-containment Lid (optional)					

For bio-containment add **"-BC"** to package no. above

## MIKRO 185 CUSTOM PACKAGE

	<b>Supports over 29 application configurations, please inquire for more information</b>					
	<b>Tube types include:</b> Microliter (0.2 - 0.8 mL adapter combinations are available, inquire for details), Micro Spin Columns, Pediatric					



[hettweb.com/mikro185](http://hettweb.com/mikro185)

DISTRIBUTED BY:



**Helmer Scientific**  
 14400 Bergen Boulevard, Noblesville, IN 46060 USA  
**Toll Free:** (800) 743-5637 | **Email:** sales@helmerinc.com  
[www.helmerinc.com](http://www.helmerinc.com)



Discover more solutions online at [hettweb.com](http://hettweb.com) **Toll Free:** 1(866) 370-4388 **Email:** info@hettweb.com  
 Copyright © 2016 Hettich Instruments, LP. All rights reserved. Lit. No. PS-MIKRO-185-E Rev A 01-2016.



LAB TECHNOLOGY