



CENTRIFUGES

ROTINA 380

BLOOD / CELL CULTURE / PLATES / BOTTLES

Small Body, Big Heart

High-performance benchtop centrifuge supporting a wide array of accessories for multiple applications in just one rotor. Its size and performance make it the ideal benchtop centrifuge for most clinical, research and industrial laboratories.

MAX. CAPACITY	4 x 290 mL
MAX. RPM/RCF	15,000/24,400
TEMP. CONTROL	from -20 to +40°C (380 R)

Distributed by:



Model ROTINA 380 without temperature control pictured.



QUIET OPERATION

Quiet as a normal conversation. Will not increase noise levels or contribute to distractions.



FIVE YEAR WARRANTY

All Hettich centrifuges are protected by our 5-year warranty and full technical support.



FUNCTIONAL DESIGN

Easy to lift lid with front-facing control panel for effortless interaction and programming.



FULLY CERTIFIED

Engineered and manufactured in accordance with all relevant safety and production standards.



ENHANCED SAFETY

Self-closing, auto-locking lid with steel latches. All-metal chamber provides reinforced security.



MADE IN GERMANY

German engineered and manufactured to be quiet, reliable and safe since 1904.

DIMS. (HxWxD): 16.5" x 18" x 23.6" (380) | 16.5" x 18" x 29.5" (380 R)

→ hettweb.com/rotina380

HETTICH ROTINA 380 APPLICATION PACKAGES



ROTINA 380 BLOOD TUBE PACKAGE - 4

	Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF	PKG. NO.	PKG. NO. (TEMP. CONTROL)
	64	1.6 - 7	13 x 75/100	≤ 5,000	≤ 4,668	380BLOOD4	380RBLOOD4
	52	4 - 10	16 x 75/100	≤ 5,000	≤ 4,668		
Includes: (1) Centrifuge, (1) Horizontal Rotor, (4) Buckets, (4) Adapters 13 x 75/100 mm, (4) Adapters 16 x 75/100, (4) Bio-containment Lids (optional)							
For bio-containment add "-BC" to package no. above							

ROTINA 380 BLOOD TUBE PACKAGE - 2

	Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF	PKG. NO.	PKG. NO. (TEMP. CONTROL)
	32	1.6 - 7	13 x 75/100	≤ 5,000	≤ 4,668	380BLOOD2	380RBLOOD2
	26	4 - 10	16 x 75/100	≤ 5,000	≤ 4,668		
Includes: (1) Centrifuge, (1) Horizontal Rotor, (4) Buckets, (2) Adapters 13 x 75/100 mm, (2) Adapters 16 x 75/100, (4) Bio-containment Lids (optional)							
For bio-containment add "-BC" to package no. above							

ROTINA 380 CELL CULTURE PACKAGE - 4

	Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF	PKG. NO.	PKG. NO. (TEMP. CONTROL)
	36	15	17 x 120	≤ 5,000	≤ 4,863	380CELLCULTURE4	380RECELLCULTURE4
	16	50	29 x 115	≤ 5,000	≤ 4,863		
Includes: (1) Centrifuge, (1) Horizontal Rotor, (4) Buckets, (4) Adapters 15 mL, (4) Adapters 50 mL, (4) Bio-containment Lids (optional)							
For bio-containment add "-BC" to package no. above							

ROTINA 380 CELL CULTURE PACKAGE - 2

	Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF	PKG. NO.	PKG. NO. (TEMP. CONTROL)
	18	15	17 x 120	≤ 5,000	≤ 4,863	380CELLCULTURE2	380RECELLCULTURE2
	8	50	29 x 115	≤ 5,000	≤ 4,863		
Includes: (1) Centrifuge, (1) Horizontal Rotor, (4) Buckets, (2) Adapters 15 mL, (2) Adapters 50 mL, (4) Bio-containment Lids (optional)							
For bio-containment add "-BC" to package no. above							

ROTINA 380 PLATE PACKAGE

	Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF	PKG. NO.	PKG. NO. (TEMP. CONTROL)
	10	Microtiter	--	≤ 4,000	≤ 2,397	380PLATE	380RPLATE
Includes: (1) Centrifuge, (1) Horizontal Rotor, (2) Buckets							

ROTINA 380 BOTTLE PACKAGE

	Bottle Qty.	Volume (mL)	Size (mm)	RPM	RCF	PKG. NO.	PKG. NO. (TEMP. CONTROL)
	4	250	61 x 122	≤ 5,000	≤ 4,863	380BOTTLE	380RBOTTLE
Includes: (1) Centrifuge, (1) Horizontal Rotor, (4) Buckets, (4) Adapters, (4) Bio-containment Lids (optional)							
For bio-containment add "-BC" to package no. above							

ROTINA 380 CUSTOM PACKAGE

	Supports over 234 application configurations, please inquire for more information						
	Tube types include: Blood, Pediatric, Urine, Conical and Skirted up to 50 mL, Standard up to 250 mL, Plates, Bottles up to 290 mL						



hettweb.com/rotina380

DISTRIBUTED BY:



Helmer Scientific
 14400 Bergen Boulevard, Noblesville, IN 46060 USA
Toll Free: (800) 743-5637 | **Email:** sales@helmerinc.com
www.helmerinc.com



Discover more solutions online at hettweb.com **Toll Free:** 1(866) 370-4388 **Email:** info@hettweb.com
 © 2016 Hettich Instruments, LP. All rights reserved. Lit. No. PS-ROTINA-380-E Rev A 01-2016.



LAB TECHNOLOGY