

ROTOFIX 32 A

BLOOD / CELL CULTURE / URINE / CYTOLOGY

Rugged and Indispensable

The ROTOFIX 32 A is a rugged, versatile centrifuge for routine laboratory tasks in doctors' practices and small hospitals. Equipped for several clinical chemistry applications including cytology. With a wide variety of accessories the ROTOFIX 32 A is also used for the preparation of samples in industrial and research laboratories.

MAX. CAPACITY	4 x 100 mL / 6 x 85 mL
MAX. RPM/RCF	6,000/4,226

Distributed by:



QUIET OPERATION

Quiet as a normal conversation. Will not increase noise levels or contribute to distractions.



FIVE YEAR WARRANTY

All Hettich centrifuges are protected by our 5-year warranty and full technical support.



FUNCTIONAL DESIGN

Easy to lift lid with front-facing control panel for effortless interaction and programming.



FULLY CERTIFIED

Engineered and manufactured in accordance with all relevant safety and production standards.



ENHANCED SAFETY

Self-closing, auto-locking lid with steel latches. All-metal chamber provides reinforced security.




MADE IN GERMANY

German engineered and manufactured to be quiet, reliable and safe since 1904.


HETTICH ROTOFIX 32 A APPLICATION PACKAGES




ROTOFIX 32 A BLOOD TUBE PACKAGE

	Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF	PKG. NO.
	20	1.6 - 7	13 x 75/100	≤ 4,000	≤ 2,558	32ABLOOD
	16	4 - 10	16 x 75/100	≤ 4,000	≤ 2,522	
Includes: (1) Centrifuge, (1) Horizontal Rotor, (4) Buckets, (4) Adapters 13 x 75/100 mm, (4) Adapters 16 x 75/100 mm, (4) Bio-containment Lids (optional)						For bio-containment add “-BC” to package no. above


ROTOFIX 32 A URINALYSIS PACKAGE

	Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF	PKG. NO.
	12	12 KOVA	17 x 102	≤ 4,000	≤ 2,451	32AURINALYSIS
	Includes: (1) Centrifuge, (1) Horizontal Rotor, (4) Buckets, (4) Adapters KOVA					


ROTOFIX 32 A CELL CULTURE PACKAGE - 4

	Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF	PKG. NO.
	16	15 Conical	17 x 120	≤ 4,000	≤ 2,612	32ACELLCULTURE4
	4	50 Conical	29 x 115	≤ 4,000	≤ 2,576	
Includes: (1) Centrifuge, (1) Horizontal Rotor, (4) Buckets, (4) Adapters 15 mL, (4) Adapters 50 mL						


ROTOFIX 32 A CELL CULTURE PACKAGE - 6

	Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF	PKG. NO.
	6	15 Conical	17 x 120	≤ 4,000	≤ 2,701	32ACELLCULTURE6
	6	50 Conical	29 x 115	≤ 4,000	≤ 2,701	
Includes: (1) Centrifuge, (1) Horizontal Rotor, (6) Buckets, (6) Adapters 15 mL						


ROTOFIX 32 A CYTOLOGY PACKAGE - 12

	Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF	PKG. NO.
	12	0.4 - 6	--	≤ 2,000	≤ 438	32ACYTOCLIP12
	Includes: (1) Centrifuge, (1) Basic Cytology Rotor, (2) Cyto Clip 6-packs, (1) Cell Funnel, disposable 50-pack					

ROTOFIX 32 A CYTOLOGY PACKAGE - 12-S

	Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF	PKG. NO.
	12	0.4 - 6	--	≤ 2,000	≤ 438	32ACYTOCLIP12-S
	Includes: (1) Centrifuge, (1) Shielded Cytology Rotor, (2) Cyto Clip 6-packs, (1) Cell Funnel, disposable 50-pack					

ROTOFIX 32 A CUSTOM PACKAGE

	Supports over 210 application configurations, please inquire for more information					
	Tube types include: Microliter up to 5 mL, Blood, Pediatric, Urine, Conical and Skirted up to 50 mL, Standard up to 100 mL, Cytology					



hettweb.com/rotofix32a

DISTRIBUTED BY:



Helmer Scientific
14400 Bergen Boulevard, Noblesville, IN 46060 USA
Toll Free: (800) 743-5637 | **Email:** sales@helmerinc.com
www.helmerinc.com



Discover more solutions online at hettweb.com **Toll Free:** 1(866) 370-4388 **Email:** info@hettweb.com
© 2016 Hettich Instruments, LP. All rights reserved. Lit. No. PS-ROTOFIX-32-A-E Rev A 01-2016.

Hettich
LAB TECHNOLOGY