

BD Pyxis™ ES Refrigerator

Secure, single-line item storage of high-risk and high-value refrigerated medication

What fits in a standard/half-height bin

| Medication type* | Quantity* |
|---|-----------|
| Vials 1 mL | 175 |
| Vials 2 mL | 100 |
| Vials 4 mL | 75 |
| Vials 5 mL | 70 |
| Vials 10 mL | 24 |
| Syringes 1 mL | 50 |
| Syringes 2 mL | 30 |
| Syringes 4 mL | 20 |
| Syringes 4 mL plus needles in separate pkg. | 10 |
| Syringes 10 mL | 12 |
| Ampules 1 mL | 40 |
| Suppositories | 75 |
| Eye drops 15 mL bottles | 30 |
| Rabies vaccine kits | 4 |

*Medication packaging may vary by manufacturer. Represented quantity is an approximation.



Standard/half-height bin internal dimensions (*w x h x d*) 4.85 in. x 1.52 in. x 12.57 in.



What fits in a deep/full-height bin

| Medication type* | Quantity* |
|--|-----------|
| Syringes 10 mL | 40 |
| PCA syringes 30 mL | 16 |
| PCA syringes 40/50 mL | 10 |
| IV bags 50 mL (with overwrap packaging) | 6 |
| IV bags 50 mL (without overwrap packaging) | 20 |
| IV bags 100 mL (with packaging) | 4 |
| IV bags 100 mL (without packaging) | 10 |
| IV bags 250 mL | 8 |
| IV bags 500 mL | 4 |
| IV bags 1 liter | 2 |
| Bottles 250/300 mL | 4 |
| Bottles 500 mL | 2 |

*Medication packaging may vary by manufacturer. Represented quantity is an approximation.



Deep/full-height bin internal dimensions (w x h x d) 4.85 in. x 4.52 in. x 12.57 in.

Common types of refrigerated medications

- Insulin
- Vaccines (e.g., hepatitis, influenza, rabies)
- Antibiotics (e.g., penicillin)
- IV solutions
- Suppositories
- Neuromuscular blockers and paralytics (e.g., Rocuronium/Zemuron™, succinylcholine chloride/Quelicin®)
- Lorazepam/Ativan®
- Compounded IV solutions
- Oncology drugs (e.g., octreotide/Sandostatin®, TBO-filgrastim/Granix™)
- Famotidine/Pepcid®
- Alteplase/Cathflo® Activase®
- Darbepoetin alfa/Aranesp®
- Oral antibiotic suspensions (e.g., amoxicillin clavulanate/Augmentin®, cephalexin/Keflex®)

BD, San Diego, CA, 92130, U.S.

bd.com

BD, the BD Logo and Pyxis are trademarks of Becton, Dickinson and Company or its affiliates. All other trademarks are the property of their respective owners. © 2019 BD. All rights reserved. 191-1691 (0319/PDF) BD-9725.

