

Installation & Operational Qualification Services

from Helmer Scientific

Validation and Regulatory Compliance

Helmer recognizes the need to achieve validation to help maintain regulatory compliance. Validating equipment through Installation and Operational Qualification ensures reliable operational criteria are met.

We offer comprehensive Installation and Operational Qualification services. Following the completion of services, thorough documentation will be available for your records.



What is Installation Qualification?

Helmer offers comprehensive Installation Qualification services that provide a documented verification process ensuring the equipment has been properly delivered, installed and configured. This offers validation that the equipment can operate as intended to standards set by the manufacturer.

IQ requirements include:

- Checking for proper location
- Checking for proper power supply
- Checking for acceptable environmental conditions
- Verifying software installation
- Verifying connections
- Recording calibrations and validation dates

What is Operational Qualification?

Operational Qualification provides verification that the equipment operate as specified by the manufacturer. The equipment is tested according to manufacturer operational guidelines.

OQ should be performed:

- After installation
- In accordance with standard operating procedures provided by the manufacturer

Contact Technical Service!

800.743.5637 | +1-317-773-9073
email: scoordinator@helmerinc.com