

Door Handles and Adapters for Pyxis® Locks

Helmer Scientific MedLock Solutions

Special accessories are available for facilities that want to install the Pyxis Remote Manager or Smart Remote Manager on Helmer upright and undercounter refrigerators.

Pyxis Compatible Door Handles for Upright Models

The door handle features a cutout and spacer which enable the Pyxis lock to wrap around the front and side of the door.

- Available for single door upright refrigerators 13, 20 and 25 cu ft (377, 572 and 714 L)
- Compatible with glass and solid doors
- Factory installed or Field Installation Kits available


Pyxis Lock Adapter Kit for Undercounter Models

The adapter kit includes an extension bar and alternate handle that mount to the door of the unit, allowing the Pyxis lock to be properly installed.

- Adapter kits available for undercounter 5 cu ft (150 L) refrigerators
- Compatible with glass and solid doors, powder-coated and stainless steel exteriors
- Factory installed or Field Installation Kits available



Cut-out of Door Handle for Upright Refrigerators



Upright Pyxis®
Compatible
Door Handle*



Undercounter Adapter*
(Shown with Pyxis® Lock)

*Confirm dimensions of your Pyxis lock prior to placing order

GX Solutions i.Series® & Horizon Series™ Models

Factory Installed Options

Pyxis Lock Adapter Kit, Factory installed, 5 cu ft (150 L)

Pyxis Lock Compatible Door Handle, Factory installed, 13 cu ft (377 L)

Pyxis Lock Compatible Door Handle, Factory installed,
20, 25 cu ft (572, 714 L)

Field Installation Kits

Pyxis Lock Adapter Kit, Field installation, 5 cu ft (150 L)

Pyxis Lock Compatible Door Handle, Field installation, 13 cu ft (377 L)

Pyxis Lock Compatible Door Handle, Field installation,
20, 25 cu ft (572, 714 L)

Note:

Adapter kits and compatible door handles do not include the medication dispensing lock.

Pyxis compatible door handles and lock adapter kits also available for legacy refrigerator models. Contact your sales rep for more information

Adjustable, hasp available from BD (BD part number 133480-01) for use on Helmer GX Solutions undercounter freezer.