


Bridge Shelf

for Double Door Refrigerators

Increase Your Usable Shelf Space


A bridge shelf spans the gap between shelves in a double-door refrigerator. With bridge shelves installed, a Helmer double-door refrigerator offers over 13% more storage area per shelf.

- Available for 45 cu ft and 56 cu ft models
- Attaches to existing shelves (no additional hardware required)
- Bacteria-resistant powder coated, galvanized steel



Increased Shelf Space

Model Size	Length	Width
45 cf GX Unit	11 in	17.3 in
45 cf GX Unit	27.94 cm	43.94 cm
56 cf GX Unit	11 in	23.3 in
56 cf GX Unit	27.94 cm	59.18 cm



GX Solutions

Bridge Shelf, 45 cu ft refrigerator, 11" w x 17.3" d
Part # 323081-1-069

Bridge Shelf, 56 cu ft refrigerator, 11" w x 23.3" d
Part # 323103-1-069

Legacy Models

Bridge Shelf, 45 cu ft refrigerator, 22" w x 18" d
Part # 400845-1

Bridge Shelf, 56 cu ft refrigerator, 22" W x 24" d
Part # 400845-2