

# Configuration Guide

## for GX Solutions Refrigerators

This configuration guide provides information on standard and maximum configurations and maximum locking drawers for GX Pharmacy, Blood Bank and Laboratory Refrigerators.

### GX Solutions - Pharmacy Models

Model	Capacity cu ft   L	Standard Configuration	Maximum Configuration*	Locking Drawer Max
iPR105-GX / HPR105-GX	5.3   150	2 ventilated drawers	2 ventilated drawers	2
iPR113-GX / HPR113-GX	13.3   377	1 ventilated shelf 5 ventilated drawers	1 ventilated shelf 5 ventilated drawers	5
iPR120-GX / HPR120-GX	20.2   572	1 ventilated shelf 6 ventilated drawers	8 ventilated drawers	1, bottom only
iPR125-GX / HPR125-GX	25.2   714	1 ventilated shelf 6 ventilated drawers	8 ventilated drawers	1, bottom only
iPR245-GX / HPR245-GX	44.9   1271	2 ventilated shelves 12 ventilated drawers	16 ventilated drawers	2, bottom only
iPR256-GX / HPR256-GX	56   1586	2 ventilated shelves 12 ventilated drawers	16 ventilated drawers	2, bottom only



### GX Solutions - Blood Bank Models

Model	Capacity cu ft   L	Standard Configuration	Maximum Configuration*	Locking Drawer Max
iBR105-GX / HBR105-GX	5.3   150	2 solid drawers	2 drawers	2
iBR113-GX / HBR113-GX	13.3   377	5 solid drawers	5 solid drawers	5
iBR120-GX / HBR120-GX	20.2   572	7 solid drawers	8 solid drawers	1, bottom only
iBR125-GX / HBR125-GX	25.2   714	7 solid drawers	8 solid drawers	1, bottom only
iBR245-GX / HBR245-GX	44.9   1271	14 solid drawers	16 solid drawers	2, bottom only
iBR256-GX / HBR256-GX	56   1586	14 solid drawers	16 solid drawers	2, bottom only



### GX Solutions - Laboratory Models

Model	Capacity cu ft   L	Standard Configuration	Maximum Configuration*	Locking Drawer Max
iLR105-GX / HLR105-GX	5.3   150	2 ventilated shelves	3 ventilated shelves	2
iLR113-GX / HLR113-GX	13.3   377	4 ventilated shelves	6 ventilated shelves	5
iLR120-GX / HLR120-GX	20.2   572	4 ventilated shelves	8 ventilated shelves	1, bottom only
iLR125-GX / HLR125-GX	25.2   714	4 ventilated shelves	8 ventilated shelves	1, bottom only
iLR245-GX / HLR245-GX	44.9   1271	8 ventilated shelves	16 ventilated shelves	2, bottom only
iLR256-GX / HLR256-GX	56   1586	8 ventilated shelves	16 ventilated shelves	2, bottom only



\*Other configurations available, not to exceed maximum configuration. Contact Helmer Scientific.