

# Pedestals

For GX Solutions Undercounter Refrigerators and Freezers

Pedestals offer greater versatility for your undercounter refrigerator or freezer. They are designed for raising units to a more comfortable working height. Medical pedestals are for use in hospitals, clinics, private health practices, pharmacies, and more.



## Undercounter Pedestal

Base elevates undercounter refrigerators and freezers above the floor to a more comfortable working height. Brackets secure models to pedestal for secure installation. Feet with plastic protectors preserve surfaces.

May be used with:

Refrigerator models iLR105-GX, HLR105-GX, iPR105-GX, HPR105-GX, iBR105-GX, HBR105-GX

Freezer models iLF105-GX, HLF105-GX, iBF105-GX, HBF105-GX

Pedestal (w x d x h):

24 x 24 x 17.5 in 610 x 610 x 445 mm

Weight: 42 lb 19 kg

Pedestal height with GX undercounter refrigerator or freezer with leveling feet\*: 48.9 in 1242mm

\*Add 1.6 in (40.7 mm) if leveling feet are fully extended

