

Locking Drawers

Locking drawers offer isolated storage for prescription and non-prescription medications and pharmaceuticals, vaccines, blood products, and more for patient care in a hospital, clinic, nursing facility, or other similar health care facility.

A ventilated cover is secured with a key lock, adding an extra layer of protection while maintaining excellent temperature uniformity for the safe storage of products.

Locking drawers are available for undercounter and upright refrigerators.



iPR256-GX shown with 2 locking drawers

Locking Drawers for GX Solutions Refrigerators

- Available for i.Series® and Horizon Series™ refrigerator models (5-56 cu ft)
- Includes key lock for added security
- Solid stainless steel pull-out drawer includes ventilated cover for optimized temperature performance
- Scratch-guard edge protectors help prevent scratches on glass doors



Specifications

	Drawer Dimensions w x h x d / in mm	Drawers per unit* place in bottom of unit
5 cu ft	16.7 x 4 x 19.2 424 x 104 x 104	2
13 cu ft	16.7 x 4 x 19.2 424 x 104 x 104	5
20 cu ft	21.5 x 4 x 20.7 546 x 104 x 527	1 / bottom position
25 cu ft	21.5 x 4 x 26.7 546 x 104 x 679	1 / bottom position
45 cu ft	21.5 x 4 x 20.7 546 x 104 x 527	2 / one per side in bottom position
56 cu ft	21.5 x 4 x 26.7 546 x 104 x 679	2 / one per side in bottom position

*Additional drawers and configurations may be possible. Please contact Helmer for more information on locking drawer options.