

Pharmacy / Vaccine Product Guide



HOSPITAL PHARMACY
HEALTH DEPARTMENT
OUTPATIENT CARE
PATIENT CARE AREAS
RETAIL PHARMACY



Pharmacy/Vaccine Product Portfolio



Maintain the highest quality of patient care and safety with temperature controlled storage from the name you trust.

Helmer Scientific offers many innovative products for hospital pharmacies, clean rooms, patient care and outpatient facilities. Helmer refrigerators and freezers are designed and manufactured specifically for the safe storage of pharmaceutical products. Safeguard critical and expensive medications, vaccines, reagents, samples and more with Helmer's true medical-grade performance.

Backed with industry leading warranties, Helmer pharmacy solutions meet your most demanding applications. Take the worry out of your cold storage and focus on providing high quality patient care.

Understanding Helmer Cold Storage Model Numbers

The first three letters of a Helmer part number represent the series, segment type and product type. These are followed by three numbers which represent the number of doors followed by the capacity of the unit. Any additional letters will be included after a hyphen and represent a special product type.

Model Nomenclature Guide Example iPR125-GX

Series	Type	Product	Number of doors	Capacity	Line
i.Series	Pharmacy	Refrigerator			
i	P	R	1	25	GX

GX Solutions Professional Medical-grade Pharmacy Refrigerators







i.Series® · Horizon Series™

+2°C to +10°C Factory set to +5°C

- » Powered by OptiCool™ cooling technology with variable capacity compressor (VCC) and natural hydrocarbon refrigerants
- » Maximum temperature uniformity and quick recovery
- » Ventilated pull-out drawers and shelves
- » Bacteria-resistant powder coating
- » Directional LED lighting (standard on Upright)
- » Self-closing doors with magnetic closure
- » Optional solid ballast simulator
- » Custom storage: drawers, ventilated drawers & ventilated shelves
- » Sustainable U.S. EPA and SNAP compliant natural refrigerants
- » NSF/ANSI 456 vaccine standard certified and supports CDC guidelines for vaccine storage
- » ADA compliant full-length handle; compatible with ADA compliant countertops



GX Solutions Pharmacy Refrigerator Selection Guide

Model Number	Capacity cu ft L	Series	Style	Electrical Frequency V (Hz)	Door(s)	Ventilated Drawers	Ventilated Shelves	Sound dB	Exterior Dimensions W x H x D (in/mm)	Net Weight lb kg
 iPR105-GX HPR105-GX	5.3 150	i.Series Horizon Series	Undercounter	115 (60) 220-240 (50/60)	1, Solid	2	-	42	24 x 31.5 x 26.7 610 x 796 x 678	182 83
 iPR113-GX HPR113-GX	13.3 377	i.Series Horizon Series	Upright	115 (60) 220-240 (50/60)	1, Dual-pane glass	5	1	39	24.8 x 70.5 x 26.2 628 x 1790 x 664 24.8 x 70.5 x 26.1 612 x 1790 x 663	338 154
 iPR120-GX HPR120-GX	20.2 572	i.Series Horizon Series	Upright	115 (60) 220-240 (50/60)	1, Dual-pane glass	6	1	49	29 x 79.6 x 28.2 737 x 2021 x 715 29 x 78.9 x 28.1 737 x 2003 x 712	438 199 487 221
 iPR125-GX HPR125-GX	25.2 714	i.Series Horizon Series	Upright	115 (60) 220-240 (50/60)	1, Dual-pane glass	6	1	49	29 x 79.6 x 34.2 737 x 2021 x 867 29 x 78.9 x 34.1 737 x 2003 x 865	517 235 514 234
 iPR245-GX HPR245-GX	44.9 1271	i.Series Horizon Series	Upright	115 (60) 220-240 (50/60)	2, Dual-pane glass	12	2	52	59 x 79.6 x 28 1499 x 2021 x 712 59 x 78.9 x 28.1 1499 x 2003 x 712	775 352 774 351
 iPR256-GX HPR256-GX	56 1586	i.Series Horizon Series	Upright	115 (60) 220-240 (50/60)	2, Dual-pane glass	12	2	52	59 x 79.6 x 34.1 1499 x 2021 x 865 59 x 78.9 x 34.1 1499 x 2003 x 865	827 376 817 371

GX Solutions Professional Medical-grade Laboratory Freezers (-30°C)




i.Series® · Horizon Series™

Upright -15°C to -35°C Factory set to -30°C
 Undercounter -15°C to -30°C Factory set to -30°C

- » Powered by OptiCool™ cooling technology with variable capacity compressor (VCC) and natural hydrocarbon refrigerants
- » Maximum temperature uniformity and quick recovery
- » Ventilated shelves
- » Bacteria-resistant powder coating
- » Self-closing doors with magnetic closure
- » Auto-defrost with programmable defrost cycles
- » Heated door frame
- » Custom storage: drawers, ventilated drawers, and ventilated shelves
- » Sustainable U.S. EPA and SNAP compliant natural refrigerants
- » Meets NSF/ANSI 456 vaccine storage standard requirements with NSF/ANSI 456 certification option
- » Compatible with ADA compliant countertops (undercounter)



GX Solutions Laboratory Freezer Selection Guide

Model Number	Capacity cu ft L	Series	Style	Electrical Frequency V (Hz)	Door(s)	Ventilated Shelves	Sound dB	Exterior Dimensions W x H x D (in/mm)	Net Weight lb kg
 iLF105-GX HLF105-GX	5.3 150	i.Series Horizon Series	Undercounter	115 (60) 220-240 (50/60)	1, Solid	2	46	24 x 31.5 x 26.7 610 x 798 x 678	179 82
 iLF120-GX HLF120-GX	20.2 572	i.Series Horizon Series	Upright	115 (60) 220-240 (50/60)	1, Solid	4	52	29 x 80.3 x 28.6 737 x 2039 x 725 29 x 80.3 x 28.4 737 x 2039 x 720	400 182 397 181
 iLF125-GX HLF125-GX	25.2 714	i.Series Horizon Series	Upright	115 (60) 220-240 (50/60)	1, Solid	4	52	29 x 80.3 x 34.4 737 x 2039 x 874	444 202 441 201

GX Solutions Professional Medical-grade Pharmacy Pass-Thru Refrigerators

i.Series® · Horizon Series™

+2°C to +10°C Factory set to +5°C

- » Powered by OptiCool™ cooling technology with variable capacity compressor (VCC) and natural hydrocarbon refrigerants
- » Maximum temperature uniformity and quick recovery
- » Helps meet USP<797> clean room requirements
- » Solid, stainless steel two-way, liquid-tight drawers and ventilated shelves
- » Bacteria-resistant powder coating
- » Directional LED lighting
- » Self-closing doors with magnetic closure
- » Optional solid ballast simulator
- » Custom storage options: drawers, ventilated shelves
- » Interlock doors (Optional on i.Series® models)
- » Sustainable U.S. EPA and SNAP compliant natural refrigerants
- » ADA compliant full-length handle



GX Solutions Pass-Thru Refrigerator Selection Guide

Model Number	Capacity cu ft L	Series	Style	Electrical Frequency V (Hz)	Door(s)	Solid Drawers	Shelves	Exterior Dimensions W x H x D (in/mm)	Net Weight lb kg
iPR226-GX HPR226-GX	26.4 747	i.Series Horizon Series	Pass-Thru	115 (60) 220-240 (50/60)	2, Dual-pane glass	3	3	29 x 79.6 x 35.8 737 x 2021 x 1014 29 x 79.3 x 35.8 737 x 2013 x 1014	526 239 523 237
iPR458-GX HPR458-GX	58.5 1656	i.Series Horizon Series	Pass-Thru	115 (60) 220-240 (50/60)	4, Dual-pane glass	6	6	59 x 79.6 x 35.8 1499 x 2021 x 1014 59 x 79.3 x 35.8 1499 x 2013 x 1014	868 394 862 391

Ultra-Low Temperature Freezers (-80°C)

i.Series®

-50°C to -86°C Factory set to -80°C

- » Heat Barrier System™ reduces frost
- » Insulated and structurally reinforced inner doors
- » Ice-resistant sealing surface
- » Robust, single-handed door handle
- » Optimized oil management
- » Integrated access control
- » Diagnostic data provides faster troubleshooting



-80° C Freezer Selection Guide

Model Number	Capacity cu ft L	Series	Style	Exterior Dimensions W x H x D (in mm)	Electrical Frequency V (Hz)	No of Shelves	Door(s)	Panels	Net Weight lb kg
iUF118	18 510	i.Series	Upright	32.5 x 78.2 x 37.5 826 x 1986 x 953	208/230 (60)	4	1	Vacuum insulated	622 282
iUF126	26 736	i.Series	Upright	43.8 x 78.2 x 37.5 1113 x 1986 x 953	208/230 (60)	4	1	Vacuum insulated	725 328

OPTION

Integrated Electronic Access Control

i.Series®

Ultra-Low Freezers - Standard
Refrigerators, Freezers - Optional

Integrated electromagnetic access control offers secure access to the refrigerator or freezer via the i.C³ Information Center. An electromagnetic lock is integrated into the unit and controlled by PIN entry.



- » Up to 100 unique User ID's can be set up directly from the i.C³ and programmed with an individual access PIN
- » The Access Control display can be used as an alternative Home Screen
- » The Access Control log provides an audit trail of door entries and can be downloaded
- » Key override is included

Horizon Series™

Refrigerators and Freezers - Optional

The Horizon Access Control keypad provides fail-secure digital keypad access. A lock is integrated into the unit and unlocked by entering a user code on the keypad



- » Up to 100 User Codes can be entered into the keypad
- » Key override included



Helmer Scientific

14400 Bergen Boulevard | Noblesville, IN 46060 USA
Toll Free (U.S. and Canada): 800.743.5637
Phone: 317.773.9073 | Fax: 317.773.9082
www.helmerinc.com

