

Protective Shelf Guards and Dividers

for Platelet Storage Systems



Helmer Scientific Solutions

Platelets are a critical blood therapy and Helmer Scientific platelet storage systems are designed to safeguard these products.

Some platelet bags utilize large amounts of tubing, are oversized, or include multiple satellite bags. For these situations, some customers have requested guards to help ensure these types of platelet bags stay in place during agitation.

Helmer has created a drawer edge and divider guard accessory to help provide additional support to secure these platelet bags.

Platelet Edge Guard & Divider Kits

Designed to prevent bag damage that could be caused from bag sets shifting or being improperly loaded on shelves.


- Guards include a magnetic tape strip which adheres to the metal drawers providing protection along the sides and back of the drawer.
- Additional kits can be used to separate storage spaces within the drawer.

Multiple Kits Available

- Compatible kits vary by model.
- Kits available for agitator model PF96-Pro, as well as for legacy models PF96i and PF96h.
- Kits available for incubator model PC4200-Pro, as well as for legacy models PC4200i and PC4200h.



Side Guard Kit

Left and right side guards designed to ensure platelets, tubing or satellite bags do not interfere with side rails of agitator drawer. Use multiple kits to partition sections of the shelf




Frame Guard Kit

This kit includes left and right side guards and front and back guards to create a guard frame around the entire shelf.



Pro Line Platelet Agitator drawer shown with Side Guard Kit



Pro Line Platelet Agitator drawer shown with Frame Guard Kit



Pro Line Platelet Agitator drawer shown with added guards

Contact the Helmer Custom Products Division to inquire about additional size availability

380475-1/B 11-23 © 2020-23 Helmer Inc. All rights reserved.



www.helmerinc.com | 800.743.5637